



# 1. Cable connection




- Adjust 'Gain' to make Mn K $\alpha$  signal (5.895 keV) around 2 V on oscilloscope. (Full range : 10 V)

## 2. Set up at 8ID-E

- (1) Canberra AIM Multichannel Analyzer → EPICS MCA display
- (2) MCA with ROIs → IDL MCA display
- (3) ADC 1 →




# 2-1. ADC




Make sure of  
4096 channels

Adjust LLD to eliminate undesirably high intensities at low channels

## 2-2. EPICS MCA display




# 2-3. IDL MCA display



# 2-3 continued

- (1) File -> Foreground -> Open detector -> 8side:mca1
- (2) ROIs (left panel) -> Clear All
- (3) Set ROIs for Mn K $\alpha$  (taller peak) using mouse buttons  
mouse middle button : left marker of ROI  
right : right  
left : setting this ROI
- (4) Put a name 'Mn Ka' for this ROI
- (5) Repeat (3) and (4) for Mn K $\beta$  (shorter peak)
- (6) Control -> Calibrate energy -> compute calibration



## 2-4. Work with real signals

Repeat ROI setup for real signals

- Clear all previous ROIs

- Typically


  - (1) “elastic” for elastic x-rays

  - (2) “fluo” for fluorescence

  - (3) “full” for full mca

  - (4) see also these ROIs in [EPICS MCA display](#)

- Unfortunately, the order of ROIs is determined by the counts of each ROI (see how the order of ROIs changes in [EPICS MCA display](#))



### 3. Data acquisition in SPEC

(1) Start SPEC: `topaz% spec8IDE`

(if not working, type `/home/users1/S8SPEC/bin/spec8IDE`)

- Note that `'/home/user1/S8SPEC/macros/common/epics_mca_3.3.mac'`  
is assumed to be installed automatically.

(2) `SPEC> qdo /home/users/8IDEUSER/local_macros/mcadet_3.mac`

(3) To activate MCA counters in SPEC,

`SPEC> mcadet_use 1` ( `mcadet_use 0` to inactivate)

(4) `SPEC> ct 10`

- Note that in “`mcadet_3.mac`” ‘`mca`’ always reads the 1st ROI

in EPICS MCA panel, ‘`sca1`’ the 2nd ROI, and ‘`sca2`’ the 3rd ROI.

(5) To save a whole MCA spectrum,

`SPEC> mca_setup`

`SPEC> getandsave_mca`

- Make sure that a filename for SPEC was assigned.



## 3-1. Read MCA in C-PLOT

- To read a MCA spectra save in SPEC, use C-PLOT  
C-PLOT> f1 scans.4

.....

Normalize points ? No

Get MCA data if present ? Yes

....

Note that **x-axis** denotes **energy in keV**, which has been calibrated in IDL MCA panel.