

Micos Rotation Stage with Pollux Driver
2011.07.08 Joe Strzalka

Hardware: Micos RS-40 in-vacuum rotation stage

Micos Pollux Controller: SMC Pollux

Adapter: Vacuum DB-9 to DB-9 for Pollux (could be eliminated in future with vacuum feedthru rewired)

Power supply for Pollux


RJ45-DB-9 cable

Software:


Joe developed a software driver to interface the Pollux driver with EPICS

Joe noticed that our medm startup script, `~/local_macros/start_epics_8idguser_current` was not running the current version of epics. He changed it to start `synApps_5_4_2` instead of `synApps_5_2_1`

EPICS:



Click "Motors 2"



Currently using 8idg:m61, or slot 1-D on the upper SBS Technologies unit on back of the 8-ID-G rack.

These are the correct parameters for the motor:

motorx_all.adl (8idg:m61) EGU: degrees
 Micor Pollux

vacuum phi

Drive	User	Dial	Raw
Hi limit	100.000000	100.000000	
Readback	-16.000000	-16.000000	-16000000
MoveAbs	-16.000000	-16.000000	-16000000
Lo limit	-100.000000	-100.000000	
MoveRel	0.000000	JogR JogF	Enable
Tweak	< 1.0000 >	HomR HomF	Disable

Stop
Pause
Move
Go

Calibration Variable

Cal Use Set Off 0.000000 Dir Pos Neg

Dynamics	Normal	Backlash	Jog	EGU
Maximum Spd	5.000000			EGU/sec
Speed	4.000000	4.000000	1.000000	EGU/sec
Base Speed	0.000000			EGU/sec
Accel. sec	0.200000	0.200000	5.000000	EGU/sec ²
Backlash distance		0.010000		EGU
Move Fraction		1.000000		

Servo	Proportional	Integral	Derivative
	0.000000	0.000000	0.000000

Resolution


Motor resolution	0.000001	EGU/step
Encoder res.	0.000001	EGU/step
Readback res.	0.000000	EGU/step
Retry deadband	0.000010	EGU
Retries 0	max: 10	
Use Encoder	No Yes	
Use Readback	No Yes	
Readback Delay	0.000000	sec
RBV inLink		
Pre-move string		
Post-move string		
Mode	supervisory	

STATUS NO_ALARM

State	0x 0x3
CurrDir	1
Moving	0
At Home	0
MotorPos	-16000000
Encoder	-16000000
MIP	0x 0x0
Err	0.000000
Version	6.44
VME Card#	0
Precision	6
Torque	Disable Enable
FWD LINK	0

Scan Ld Go Abort More

V2.6



NOTE:

- 1) This rotation stage has a home switch.
- 2) The Pollux driver treats the home position as a limit switch, so it is not possible to drive the motor past the home position.
- 3) Joe set up the EPICS driver in such a way that the user and dial positions can be zeroed, but they cannot be set to arbitrary values.
- 4) When the Pollux driver powers up, it calls the current position zero.

More details: http://www.micos-online.com/web2/en/1,2,060,smc_pollux.html
<http://www.micos-online.com/web2/en/1,5,110,rs40.html>